



LEGENDScreen[™]

Lyophilized Antibody Array

Catalog Number: 700011

Size: 1 kit, ready to use
1 test per antibody/well

Reactivity: Human

Antibody Format: PE-conjugated, lyophilized

Configuration: 364 pre-titrated antibodies, including 354 specificities and 10 isotype controls, arrayed on four 96-well plates, with one specificity per well

It is highly recommended that this manual be read in its entirety before using this product.

Do not use this kit beyond the expiration date

For research use only

BioLegend, Inc
Biolegend.com

For Research Purposes Only. Not for use in diagnostic or therapeutic procedures. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of BioLegend is strictly prohibited.

Table of Contents

Description	4
Materials Provided	4
Materials to be Provided by the End User	5
Storage and Handling	5
Preparation of Cells for Staining	5
Plate Preparation	6
Cell Staining Procedure	7
Tips for Successful Staining	7
Data Analysis	8
Frequently Asked Questions	9
Product Performance	10
Antibody information Tables	15

Description

The LEGENDScreen™ Human PE Antibody kit contains 364 PE-conjugated monoclonal antibodies to cell surface markers as well as 10 mouse, rat, or hamster Ig isotype controls in a convenient 96 well, ready to use format. The kit can be used for screening human leukocytes, cell lines and cells isolated from tissues. If other co-stains are needed, they can be added to the wells following reconstitution. Positive “hits” from the screening can be quickly identified based on the plate map (Appendix I) and the catalog number and clone information can be obtained using the antibody information tables (Appendix II). Analysis with individual fluorochrome-conjugated antibodies should be performed to confirm the screening results. The LEGENDScreen™ Human PE Antibody Kit provides a convenient, easy to use, and powerful tool for immunology, stem cell, and cancer research.

Materials Provided

Description	Quantity	Size	Part#
LEGENDScreen™ Human PE Kit Plate #1	1 Plate	96-well	750002447
LEGENDScreen™ Human PE Kit Plate #2	1 Plate	96-well	750002448
LEGENDScreen™ Human PE Kit Plate #3	1 Plate	96-well	750002449
LEGENDScreen™ Human PE Kit Plate #4	1 Plate	96-well	750002450
Cell Staining Buffer	1 Bottle	500 mL	420201
Fixation Buffer	1 Bottle	100 mL	420801
Plate Sealers (Clear)	12 Sheets	79.4mm x 141mm	78101

Additional Notes on Kit Contents:

Depending on the lot received, this kit may include additional liquid antibody vials with further instructions for their use. Not all lots will include additional liquid vials. Please thoroughly check package contents for additional vials.

To verify kit lot content, please contact our technical service team (tech@biolegend.com) with the kit Lot# provided on the product box for information.

Materials to be Provided by the End User

1. Adjustable multichannel pipettes for measuring volumes ranging from 25µL to 1,000µL
2. Centrifuge with a rotor and adaptors for 96-well plates
3. Cell culture medium (for cell culture)
4. Cell dissociation buffer (for adherent cells)
5. 1X PBS (Phosphate-Buffered Saline): 8.0 g NaCl, 1.16 g Na₂HPO₄, 0.2 g KH₂PO₄, 0.2 g KCl, add deionized water to 1 liter; pH to 7.4, 0.2µm filtered
6. Plastic reservoirs for pipetting deionized water or Cell Staining or Fixation Buffer with a multichannel pipette
7. A flow cytometer, preferably compatible with reading 96-well plates

Storage and Handling

1. Store unopened kit components at 2 - 8°C. Do not use this kit beyond its expiration date.
2. Once opened, reconstitute plates with distilled, deionized water. The reconstituted plates can be used immediately or sealed for storage at 2 - 8°C in the dark for up to one month.
3. Keep the buffer(s) at 2 - 8°C and use within one month after opening.

Kit protocol:


Preparation of Cells for Staining

1. Obtain desired tissue (*e.g. leukocytes, cell lines or tissue isolates*) and prepare a single cell suspension. Wash cells in 1X PBS or cell culture medium of choice, and resuspend in Cell Staining Buffer at a density between 4×10^6 and 1×10^7 cells/mL.
2. For cultured cells in suspension, spin and resuspend cells in Cell Staining Buffer at a density between 4×10^6 and 1×10^7 cells/mL.
3. For cultured adherent cells, dissociate cells using a mild enzyme or non- enzymatic dissociation buffer. Wash cells in 1X PBS or cell culture medium of choice, and resuspend in Cell Staining Buffer at a density between 4×10^6 and 1×10^7 cells/mL.
4. Filter the cells through a 40µm cell strainer to remove any clumps. Keep the cells on ice before use. Calculate the total volume of cell suspension needed for all the wells that will be used. Each well needs 75µL of cell suspension (~ $3 - 7.5 \times 10^5$ cells/well). Plate maps can be found in Appendix I.
 - Lower cell density (*e.g. 1.5×10^6 cells/mL*) may be used depending on the application.

Plate Preparation

One hour before the staining, perform the following steps to prepare the plates:

1. Remove the lyophilized plates from the aluminum pouches. The foil seal can be easily pierced. Handle plates with care by holding along the sides of the plate, once they have been removed from the pouches.
2. Centrifuge the plates at $600 \times g$ for 10 minutes.
3. Ensure the lyophilized cakes have settled to the bottom of the plates. Keep the plates upright at all times from this point forward. Handle plates with care so that the cakes are not agitated at any time.
4. Fill a reservoir (Fisher Scientific Cat# 14-387-069 or equivalent) with deionized water. To reconstitute the lyophilized antibodies, fill the multi-channel pipette with $25\mu\text{L}$ /well of deionized water and pierce the plate seal at a 45° angle (see images below for example of good angle and bad angle for piercing the plate seal). This makes it easier to pierce the plate seal. At a 90° angle, more force is needed to pierce the seal.



5. Dispense $25\mu\text{L}$ of deionized water **very slowly** to each well, being careful not to cause well contents to splash back onto the film. To avoid cross-contamination, discard the tips after one row or column is reconstituted. Incubate the plate at ambient temperature for 10-15 minutes before removing the seal.
6. Hold the plate firmly then slowly and carefully remove the plate seal, starting from one corner, uncovering one well at a time, for easy opening and to prevent cross-well contamination. You may need to wiggle the seal a little while pulling back slowly.
7. **Discard the plate seal.**
8. Proceed with staining procedure (see next section).
9. If staining cannot be performed immediately, seal the plates with the clear plate sealers provided in the kit.
10. Keep the plates in dark until ready to use.
11. *(If necessary)*— depending on the lot, you may have received additional liquid antibody vials for verification purposes. Please follow provided instructions to add appropriate antibody amounts to blank wells on plate 4. If your kit did not come with additional liquid antibody vials, skip this step

Notes:

- PE-conjugated antibodies are light-sensitive. Try to minimize the exposure of the plates to light as much as practically possible.
- Do not open the pouches until the day you are ready to run the experiment. Once the plates are removed from the pouches, the antibodies must be reconstituted immediately.
- If an experiment is not performed after reconstitution, plates can be sealed and stored in the dark at $2 - 8^\circ\text{C}$ for up to one month.

Cell Staining Procedure

1. Using a multichannel pipette, add 75µL of cells ($\sim 3 - 7.5 \times 10^5$ cells/well) to each well of the plates.
2. Set up extra tubes to stain cells for flow cytometer setup and compensation, if needed.
3. Using the multichannel pipette, gently mix the cell suspensions by pipetting up and down 2 - 3 times. **Be sure to change tips between each row or column.** Avoid creating bubbles while pipetting.
4. Incubate for 20 - 30 minutes at 2 - 8°C in the dark.
5. Spin the plates at $500 \times g$ for 10 minutes to pellet cells. Immediately after centrifugation, dump the supernatant by quickly inverting and flicking the plate. Gently blot the plate on a clean paper towel, being careful not to disturb the cell pellet.
6. Using a multi-channel pipette add 200µL of Cell Staining Buffer to each well. Gently mix up and down to resuspend cells. **Be sure to change tips between each row or column.**
7. Repeat step 5.
8. To fix the cells, using a multichannel pipette, aliquot 100µL of Fixation Buffer into each well. Gently mix up and down to resuspend cells. **Be sure to change tips between each row or column.**
9. Incubate for 10 minutes at room temperature in the dark.
10. Repeat steps 5 - 6.
11. Repeat step 5 one more time. Resuspend cells completely in 160µL of Cell Staining Buffer per well and acquire on a flow cytometer. We recommend acquiring 70µL of sample and collecting 5,000 - 10,000 events. Users should determine the optimal number of events to be collected based on specific application they are testing. While the first plate is being acquired, store the other plates at 2 - 8°C in the dark.

Tips for Successful Staining

1. Read the entire manual carefully before the experiment.
2. Plan the experiment in advance. Designate a full day for this experiment. Do not rush any step.
3. Make sure that the flow cytometer's autosampler is well maintained and working well before the experiment. If the flow cytometer does not have an autosampler, the samples from each well of the plates should be transferred to individual FACS tubes and adjust the sample volume to approximately 300µL or more to avoid running dry. Alternatively, the samples from each well can be transferred to a Cluster Tube System (such as Corning catalog # 4410) and acquired manually. The cluster tube can then be transferred to a regular FACS tube for acquisition.
4. Make sure that enough cells have been prepared for the staining. If there are not enough cells, you may choose to divide the staining into two separate experiments.
5. Depending on the application, the cell number needed for staining can be decreased. Successful staining has been done with 1×10^5 cells/well.
6. Make sure to prepare cells for machine setup and compensation. These cells should be treated the same way as the cells for staining in the plates.
7. Handle the plates with care. Keep the plates upright at all times and be careful not to knock the plates over.
8. Protect the plates from exposure to light as much as possible.
9. Use care when handling the plates before reconstitution. The foil seal can be easily pierced. When reconstituting cakes ensure the pipette tips are at a 45° angle to the wells (**see Plate Preparation for instructions**).
10. Use extra caution when removing the plate seal after reconstitution. With the plate on a flat surface, firmly grasp the foil seal at one corner and slowly remove the seal, one well at a time. This will prevent well to well mixing of the reconstituted product (**see Plate Preparation for instructions**).
11. Some cell surface markers are sensitive to enzymatic digestion. If adherent cells are being used for staining, a mild enzyme or non-enzymatic dissociation buffer should be used when possible.
12. Make sure cells are in a single cell suspension. DNase treatment is recommended to avoid clumps caused by dead cells followed by filtration through a 40 µm cell strainer.
13. Acquire only 70µL of the 160µL total volume so that a second run can be performed if necessary.

Data Analysis

We recommend labeling acquisition files as Plate 1_A1, Plate 1_A2 and so on.

After sample acquisition is complete, transfer the files into FlowJo™, or equivalent analysis software. Using the “Create Group” function, group samples based on the isotype. This will help set gates for all samples of the same isotype. Antibodies are arranged by isotypes (isotype control followed by all antibodies of that isotype) – for example, Plate 1, A2 is Armenian hamster IgG isotype control and antibodies A3 through A4 are all of the Armenian hamster IgG isotype.

Create 10 groups based on the different isotype control antibodies as described below –

1. Armenian Hamster IgG isotype – Plate 1_A2-A4
2. Mouse IgG1, κ isotype – Plate 1_A5-H12, Plate 2_A2-H12 and Plate3 A2-C11
3. Mouse IgG2a, κ isotype – Plate 3_C12-H12 and Plate4_A2-A8
4. Mouse IgG2b, κ isotype – Plate4_A9-D12
5. Mouse IgG3, κ isotype – Plate 4_E1-E3
6. Mouse IgM, κ isotype – Plate 4_E4-E8
7. Rat IgG1, κ isotype – Plate 4_E9-E11
8. Rat IgG2a, κ isotype –Plate 4_E12-G1
9. Rat IgG2b, κ isotype –Plate 4_G2-G6
10. Rat IgM, κ isotype –Plate 4_G7-G8

Once the groups are created in FlowJo™, go to the first group and gate around the population(s) of interest and then make a histogram plot for PE. Using the isotype control sample, set a gate on the positive population and then apply it to all files in that group. Repeat this process for the remaining 9 groups. Add statistics for percentage of positive, median fluorescence intensity (MFI), event count and any other statistic of interest. Import into an excel file for further analysis using the table editor function.

If comparing control and test samples, then the histograms can be overlaid for visual comparison. If using software other than FlowJo™ please follow a similar strategy.

Frequently Asked Questions

Q: Why did or didn't my kit include additional liquid antibody vials?

A: Based on our stringent quality control procedures, certain lots of this kit may include additional liquid antibody vials. These are typically provided for additional data verification on certain markers/wells for analysis. If provided, follow additional instructions included along with the liquid antibodies in the kit. For further information, please contact our technical service team (tech@biolegend.com) with the kit lot information.

Q: What is the level of variability from one experiment to the other?

A: If the protocol is followed the variability should be minimal. The variability should be similar to single vial antibody staining.

Q: How should the kit be stored?

A: The kit should be stored at 2 - 8°C upon receipt. Once opened, the plates must be reconstituted immediately. Reconstituted plates can be used or stored at 2 - 8°C sealed in the dark and used within a month.

Q: How do I request a custom LEGENDScreen™ product with only my specificities of interest?

A: For more info, visit: biolegend.com/custom_solutions

Q: What are the guarantees regarding the lyophilized plate compared to the reconstituted plate?

A: Lyophilized product has a guaranteed shelf life of 6 months unopened. Reconstituted plates can be used or stored at 2 - 8°C sealed in the dark and used within a month. Be sure to properly seal the plates to prevent evaporation and shield the antibodies from light.

Q: I have added my own antibody solution to the lyophilized product, will the lyophilized antibody work?

A: Yes, as long as the fluorophores on these antibodies are compatible and proper compensation has been applied during acquisition and analysis.

Q: I am not going to use all the reconstituted antibody solution. Can I keep the leftover for later or re-dry the solution in the dark?

A: The antibody is in a one test per well format. There will not be any antibody left if the full test is used. Customers may decide to use less than the recommended volume per test, but this is not recommended and the performance is not guaranteed. Customers may also selectively transfer certain antibodies from the original plate to a new plate and use after reconstitution. If any antibody is not used after reconstitution, the plate can be sealed and stored at 2 - 8°C for a month in the dark. Once reconstituted, re-drying is not recommended, as this may result in a loss of signal.

Q: If I don't have enough cells and use less than 4 x 10⁶/mL (3 x 10⁵ cells/well), will it still work?

A: The kit may work with lower numbers of total cells, but we recommend trying to keep higher concentrations of cells for faster analysis. Of course, how many cells are needed depend on the specific application. Successful staining has been done with 1 x 10⁵ cells/well.

Q: Are these plates made under sterile conditions?

A: The plates are not sterile. Handle them as you would handle a typical flow cytometry staining protocol or reagent.

Q: Can I use half or less of the plate and keep the rest for later?

A: Yes. Customers can use half of the plate or whatever specificities they are interested in. However, the whole plate should be reconstituted. The half plate of antibodies must be transferred to another empty plate for the staining. The remaining half must be sealed and stored at 2 - 8°C in the dark and used within a month.


Product Performance

The LEGENDScreen™ Human PE Kit was tested and compared with BioLegend’s cataloged single vial liquid antibody reagents. For cell staining, human PBMCs were isolated and 3×10^5 cells were added to each well after the lyophilized antibodies were reconstituted. The cells were then stained for 20 minutes at 2 - 8°C, washed, and fixed with Fixation Buffer. The cells were then washed, resuspended in 160µL of Cell Staining Buffer, and analyzed using a BD FACSCanto II Flow Cytometry Analyzer.


The staining patterns and median fluorescence intensity between the lyophilized product and liquid antibodies are similar. Below are some representative data obtained when comparing the LEGENDScreen™ Human (PE) Kit vs. equivalent liquid antibodies.

Representative data

Liquid Reagents



Legendscreen Reagents



Appendix I - Plate Maps

Microsoft Excel version can be downloaded from the corresponding LEGENDScreen™ product webpage on the BioLegend website

Plate 1

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	Armenian Hamster IgG isotype control	CD278 (ICOS)	IFN- γ R β chain	Mouse IgG1, k isotype control	CD46	CD70	CD1a	CD2	β 2-microglobulin	B7-H4	Cadherin 11
B	CD10	CD100	CD103 (Integrin α E)	CD105 (Endoglin)	CD106	CD107a (LAMP-1)	CD107b (LAMP-2)	CD109	CD111 (Nectin-1)	CD112 (Nectin-2)	CD114 (G-CSFR)	CD116
C	CD119 (IFN- γ R α chain)	CD11a	CD122 (IL-2R β)	CD123	CD126 (IL-6R α)	CD127 (IL-7R α)	CD13	CD131	CD134 (OX40)	CD135 (Flt-3/Flk-2)	CD137 (4-1BB)	CD137L (4-1BB Ligand)
D	CD138 (Syndecan-1)	CD14	CD140a (PDGFR α)	CD140b (PDGFR β)	CD141 (Thrombomodulin)	CD142	CD143 (Angiotensin-converting enzyme)	CD146	CD148	CD15 (SSEA-1)	CD150 (SLAM)	CD151 (PETA-3)
E	CD154	CD156c (ADAM10)	CD158e1 (KIR3DL1, NKB1)	CD16	CD162	CD163	CD164	CD165	CD166	CD169 (Sialoadhesin, Siglec-1)	CD170 (Siglec-5)	CD172a/b (SIRP α / β)
F	CD172g (SIRP γ)	CD178 (Fas-L)	CD179a (VpreB)	CD179b (Ig λ 5)	CD18	CD180 (RP105)	CD182 (CXCR2)	CD183 (CXCR3)	CD185 (CXCR5)	CD19	CD191 (CCR1)	CD194 (CCR4)
G	CD1b	CD1c	CD200 (OX2)	CD200 R	CD202b (Tie2/Tek)	CD203c (E-NPP3)	CD205 (DEC-205)	CD206 (MMR)	CD207 (Langerin)	CD21	CD213a1 (IL-13R α 1)	CD213a2 (IL13R α 2)
H	CD218a (IL-18R α)	CD221 (IGF-1R)	CD223 (LAG-3)	CD226 (DNAM-1)	CD227 (MUC-1)	CD229 (Ly-9)	CD23	CD231 (TALLA)	CD244 (2B4)	CD245 (p220/240)	CD25	CD261 (DR4, TRAIL-R1)

Plate 2

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	CD11b	CD262 (DR5, TRAIL-R2)	DcR1 (TRAIL-R3, CD263)	CD266 (Fn14, TWEAK R)	CD268 (BAFF-R)	CD27	CD271 (NGFR)	CD275 (B7-H2, B7-RP1, ICOSL)	CD276 (B7-H3)	CD147	CD184 (CXCR4)
B	CD28	CD29	CD290 (TLR10)	CD298	CD30	CD300c	CD309 (VEGFR2)	CD31	CD314 (NKG2D)	CD317 (BST2, Tetherin)	CD324 (E-Cadherin)	CD325 (N-Cadherin)
C	CD328 (Siglec-7)	CD33	CD334 (FGFR4)	CD335 (NKp46)	CD336 (NKp44)	CD337 (NKp30)	CD34	CD340 (erbB2/HER-2)	CD344 (Frizzled-4)	CD35	CD354 (TREM-1)	CD365 (Tim-1)
D	CD366 (Tim-3)	CD367 (CLEC4A)	CD36L1 (SCARB1, SR-BI)	CD38	CD39	CD40	CD41	CD42b	CD43	CD44	CD45	CD47
E	CD48	CD49a	CD49b	CD49c (integrin α 3)	CD49d	CD5	CD50 (ICAM-3)	CD54	CD55	CD56 (NCAM)	CD58 (LFA-3)	CD6
F	CD61	CD62E	CD62L	CD62P (P-Selectin)	CD63	CD64	CD69	CD73 (Ecto-5'-nucleotidase)	CD74	CD79b (Ig β)	CD8	CD80
G	CD81 (TAPA-1)	CD82	CD83	CD85g (ILT7)	CD85k (ILT3)	CD87	CD89	CD8a	CD9	CD90 (Thy1)	CD93	CD94
H	CD95 (Fas)	CD96 (TACTILE)	CD97	CD270 (HVEM, TR2)	CXCL16	Delta-like protein 1 (DLL1)	Delta-like protein 4 (DLL4)	DR3 (TRAMP)	EGFR	GPR19	GPR56	HLA-E

Plate 3

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	Ig light chain k	IgM	Integrin $\alpha 9\beta 1$	Jagged 2	Ksp37 (FGFBP2)	LY6G6D (C6orf23)	MERTK	MSC	MSC and NPC (W4A5)	TNAP	MUC-13
B	NKp80	Notch 1	Notch 3	Notch 4	NPC (57D2)	CD352 (NTB-A)	PSMA (FOLH1)	ROR1	Siglec-10	CD328 (Siglec-7)	Siglec-8	Siglec-9
C	SSEA-5	SUSD2	α/β T Cell Receptor	TCR γ/δ	Tim-4	TLT-2	TM4SF20	TRA-2-49	TRA-2-54	TSLPR (TSLP-R)	VEGFR3 (FLT-4)	Mouse IgG2a, k isotype control
D	LOX-1	CD20	CD52	CD117 (c-kit)	CD133	APCDD1 (DRAPC1)	CD272 (BTLA)	CD198 (CCR8)	CCRL2	CD102 (ICAM-2)	CD104	CD124 (IL-4R α)
E	CD130 (gp130)	CD144 (VE-Cadherin)	CD152 (CTLA-4)	CD155 (PVR)	CD158b (KIR2DL2/L3, NKAT2)	TMEM8A	CD186 (CXCR6)	CD192 (CCR2)	CD197 (CCR7)	CD199 (CCR9)	CD209 (DC-SIGN)	CD217
F	CD230 (Prion)	CD24	CD243 (MDR-1)	CD26	CD269 (BCMA)	CD282 (TLR2)	CD284 (TLR4)	CD301 (CLEC10A)	CD303 (BDCA-2)	CD304 (Neuropilin-1)	CD307e (FcRL5)	CD323 (JAM3)
G	CD357 (GITR)	CD36	CD369 (Dectin-1/CLEC7A)	CD370 (CLEC9A/DNGR1)	CD371 (CLEC12A)	CD45RO	CD51	CD59	CD7	CD71	CD84	CD88 (CSaR)
H	Cd355 (CRTAM)	erbB3/HER-3	FPR3 (FPRL2)	Ganglioside GD2	GPR83	HLA-A,B,C	HLA-DR	Ig light chain λ	IgD	IL-28RA	integrin $\beta 5$	KLRG1 (MAFA)

Plate 4

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	MICA/MICB	SUSD2	Notch 2	TACSTD2 (TROP2)	CD3	CD99	TIGIT (VSTM3)	Mouse IgG2b, κ Isotype Control	C3AR	CCX-CKR (CCRL1)	CD11c
B	CD129 (IL-9 R)	CD158 (KIR2DL1/S1/S3/S5)	CD181 (CXCR1)	CD193 (CCR3)	CD196 (CCR6)	CD1d	CD22	CD220	CD235ab	CD258 (LIGHT)	CD274 (B7-H1, PD-L1)	CD319 (CRACC)
C	CD32	CD326 (Ep-CAM)	CD338 (ABCG2)	CD368 (CLEC4D)	CD45RA	CD45RB	CD49e	CD66a/c/e	CD85h (ILT1)	CD85j (ILT2)	CD86	CD92
D	CXCR7	Delta Opioid Receptor	Dopamine Receptor D1	EphA2	FcεRIα	GARP (LRRC32)	CD215 (IL-15Rα)	Lymphotoxin b Receptor (LT-βR)	MRGX2	CD254 (TRANCE, RANKL)	CD318 (CDCP1)	CD279 (PD-1)
E	Mouse IgG3, κ isotype control	CD255 (TWEAK)	SSEA-4	Mouse IgM, κ isotype control	CD160	CD57	CD66b	TRA-1-60-R	Rat IgG1, κ isotype control	CD115 (CSF-1R)	CD201 (EPCR)	Rat IgG2a, κ isotype control
F	CD161	CD120b	CD210 (IL-10 R)	CD267 (TACI)	CD294 (CRTH2)	CD49f	CD85a (ILT5)	CD85d (ILT4)	IgG Fc	Integrin β7	XCR1	CD360 (IL-21R)
G	Podoplanin	Rat IgG2b, κ isotype control	CD132 (common γ chain)	CD195 (CCR5)	CD4	CX3CR1	Rat IgM, κ isotype control	SSEA-3	Blank	Blank	Blank	Blank
H	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank

Appendix II - Antibody Information Tables

Microsoft Excel version can be downloaded from the corresponding LEGENDScreen™ product webpage on the BioLegend website

Plate 1

Well	Cat #	Clone	Description	Isotype
A1			Blank	
A2	400908	HTK888	Armenian Hamster IgG isotype control	Armenian Hamster IgG
A3	313508	C398.4A	CD278 (ICOS)	Armenian Hamster IgG
A4	308504	2HUB-159	IFN- γ R β chain	Hamster IgG
A5	400114	MOPC-21	Mouse IgG1, κ isotype control	Mouse IgG1, κ
A6	352402	TRA-2-10	CD46	Mouse IgG1, κ
A7	355104	113-16	CD70	Mouse IgG1, κ
A8	300106	HI149	CD1a	Mouse IgG1, κ
A9	300208	RPA-2.10	CD2	Mouse IgG1, κ
A10	395703	A17082A	β 2-microglobulin	Mouse IgG1, κ
A11	358104	MIH43	B7-H4	Mouse IgG1, κ
A12	368704	16G5	Cadherin 11	Mouse IgG1, κ
B1	312204	HI10a	CD10	Mouse IgG1, κ
B2	328408	A8	CD100	Mouse IgG1, κ
B3	350206	Ber-ACT8	CD103 (Integrin α E)	Mouse IgG1, κ
B4	800504	SN6h	CD105 (Endoglin)	Mouse IgG1, κ
B5	305806	STA	CD106	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
E1	310806	24-31	CD154	Mouse IgG1, κ
E2	352704	SHM14	CD156c (ADAM10)	Mouse IgG1, κ
E3	312708	DX9	CD158e1 (KIR3DL1, NKB1)	Mouse IgG1, κ
E4	302008	3G8	CD16	Mouse IgG1, κ
E5	328806	KPL-1	CD162	Mouse IgG1, κ
E6	333606	GHI/61	CD163	Mouse IgG1, κ
E7	324808	67D2	CD164	Mouse IgG1, κ
E8	329010	SN2 (N6-D11)	CD165	Mouse IgG1, κ
E9	343904	3A6	CD166	Mouse IgG1, κ
E10	346004	7-239	CD169 (Sialoadhesin, Siglec-1)	Mouse IgG1, κ
E11	352004	1A5	CD170 (Siglec-5)	Mouse IgG1, κ
E12	323806	SE5A5	CD172a/b (SIRP α / β)	Mouse IgG1, κ
F1	336606	LSB2.20	CD172g (SIRP γ)	Mouse IgG1, κ
F2	306407	NOK-1	CD178 (Fas-L)	Mouse IgG1, κ
F3	347404	HSL96	CD179a (VpreB)	Mouse IgG1, κ
F4	349804	HSL11	CD179b (Ig λ 5)	Mouse IgG1, κ
F5	302108	TS1/18	CD18	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
B6	328608	H4A3	CD107a (LAMP-1)	Mouse IgG1, κ
B7	354304	H4B4	CD107b (LAMP-2)	Mouse IgG1, κ
B8	323305	W7C5	CD109	Mouse IgG1, κ
B9	340404	R1.302	CD111 (Nectin-1)	Mouse IgG1, κ
B10	337410	TX31	CD112 (Nectin-2)	Mouse IgG1, κ
B11	346106	LMM741	CD114 (G-CSFR)	Mouse IgG1, κ
B12	305908	4H1	CD116	Mouse IgG1, κ
C1	308606	GIR-208	CD119 (IFN-γ R α chain)	Mouse IgG1, κ
C2	301208	HI111	CD11a	Mouse IgG1, κ
C3	339006	TU27	CD122 (IL-2Rβ)	Mouse IgG1, κ
C4	396703	S18016F	CD123	Mouse IgG1, κ
C5	352804	UV4	CD126 (IL-6Rα)	Mouse IgG1, κ
C6	351304	A019D5	CD127 (IL-7Rα)	Mouse IgG1, κ
C7	301704	WM15	CD13	Mouse IgG1, κ
C8	306104	1C1	CD131	Mouse IgG1, κ
C9	350004	Ber-ACT35 (ACT35)	CD134 (OX40)	Mouse IgG1, κ
C10	313306	BV10A4H2	CD135 (FIt-3/Flk-2)	Mouse IgG1, κ
C11	300803	S18014C	CD137 (4-1BB)	Mouse IgG1, κ
C12	311504	5F4	CD137L (4-1BB Ligand)	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
F6	312906	MHR73-11	CD180 (RP105)	Mouse IgG1, κ
F7	320706	5E8/CXCR2	CD182 (CXCR2)	Mouse IgG1, κ
F8	353706	G025H7	CD183 (CXCR3)	Mouse IgG1, κ
F9	356904	J252D4	CD185 (CXCR5)	Mouse IgG1, κ
F10	302208	HIB19	CD19	Mouse IgG1, κ
F11	362904	5F10B29	CD191 (CCR1)	Mouse IgG1, κ
F12	359412	L291H4	CD194 (CCR4)	Mouse IgG1, κ
G1	329108	SN13 (K5-1B8)	CD1b	Mouse IgG1, κ
G2	331506	L161	CD1c	Mouse IgG1, κ
G3	399803	A18042B	CD200 (OX2)	Mouse IgG1, κ
G4	329306	OX-108	CD200R	Mouse IgG1, κ
G5	334206	33.1 (Ab33)	CD202b (Tie2/Tek)	Mouse IgG1, κ
G6	324606	NP4D6	CD203c (E-NPP3)	Mouse IgG1, κ
G7	359204	HD83	CD205 (DEC-205)	Mouse IgG1, κ
G8	321106	15-2	CD206 (MMR)	Mouse IgG1, κ
G9	352204	10E2	CD207 (Langerin)	Mouse IgG1, κ
G10	354904	Bu32	CD21	Mouse IgG1, κ
G11	360404	SS12B	CD213a1 (IL-13Rα1)	Mouse IgG1, κ
G12	354404	SHM38	CD213a2 (IL13Rα2)	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
D1	356504	MI15	CD138 (Syndecan-1)	Mouse IgG1, κ
D2	399203	S18004B	CD14	Mouse IgG1, κ
D3	323506	16A1	CD140a (PDGFRα)	Mouse IgG1, κ
D4	323606	18A2	CD140b (PDGFRβ)	Mouse IgG1, κ
D5	344104	M80	CD141 (Thrombomodulin)	Mouse IgG1, κ
D6	365204	NY2	CD142	Mouse IgG1, κ
D7	344204	5-369	CD143 (Angiotensin- converting enzyme)	Mouse IgG1, κ
D8	361006	P1H12	CD146	Mouse IgG1, κ
D9	328708	A3	CD148	Mouse IgG1, κ
D10	323006	W6D3	CD15 (SSEA-1)	Mouse IgG1, κ
D11	306308	A12 (7D4)	CD150 (SLAM)	Mouse IgG1, κ
D12	350408	50-6	CD151 (PETA-3)	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
H1	313808	H44	CD218a (IL-18Rα)	Mouse IgG1, κ
H2	351806	1H7/CD221	CD221 (IGF-1R)	Mouse IgG1, κ
H3	369306	11C3C65	CD223 (LAG-3)	Mouse IgG1, κ
H4	338306	11A8	CD226 (DNAM-1)	Mouse IgG1, κ
H5	355604	16A	CD227 (MUC-1)	Mouse IgG1, λ
H6	326108	HLy-9.1.25	CD229 (Ly-9)	Mouse IgG1, κ
H7	338508	EBVCS-5	CD23	Mouse IgG1, κ
H8	329406	SN1a (M3-3D9)	CD231 (TALLA)	Mouse IgG1, κ
H9	329508	C1.7	CD244 (2B4)	Mouse IgG1, κ
H10	334403	DY12	CD245 (p220/240)	Mouse IgG1, κ
H11	356104	M-A251	CD25	Mouse IgG1, κ
H12	307206	DJR1	CD261 (DR4, TRAIL-R1)	Mouse IgG1, κ

Plate 2

Well	Cat #	Clone	Description	Isotype
A1			Blank	
A2	301306	ICRF44	CD11b	Mouse IgG1, κ
A3	307406	DJR2-4 (7-8)	CD262 (DR5, TRAIL-R2)	Mouse IgG1, κ
A4	307006	DJR3	DcR1 (TRAIL-R3, CD263)	Mouse IgG1, κ
A5	314004	ITEM-1	CD266 (Fn14, TWEAK R)	Mouse IgG1, κ
A6	316906	11C1	CD268 (BAFF-R)	Mouse IgG1, κ
A7	393211	QA17A18	CD27	Mouse IgG1, κ
A8	345106	ME20.4	CD271 (NGFR)	Mouse IgG1, κ
A9	329806	9F.8A4	CD275 (B7-H2, B7-RP1, ICOSL)	Mouse IgG1, κ
A10	351004	MIH42	CD276 (B7-H3)	Mouse IgG1, κ
A11	306212	HIM6	CD147	Mouse IgG1, κ
A12	304504	QA18A64	CD184 (CXCR4)	Mouse IgG1, κ
B1	302908	CD28.2	CD28	Mouse IgG1, κ
B2	303004	TS2/16	CD29	Mouse IgG1, κ
B3	354604	3C10C5	CD290 (TLR10)	Mouse IgG1, κ
B4	341704	LNH-94	CD298	Mouse IgG1, κ
B5	333906	BY88	CD30	Mouse IgG1, κ
B6	334804	TX45	CD300c	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
E1	336708	BJ40	CD48	Mouse IgG1, κ
E2	328304	TS2/7	CD49a	Mouse IgG1, κ
E3	359308	P1E6-C5	CD49b	Mouse IgG1, κ
E4	343803	ASC-1	CD49c (integrin α3)	Mouse IgG1, κ
E5	304304	9F10	CD49d	Mouse IgG1, κ
E6	300608	UCHT2	CD5	Mouse IgG1, κ
E7	330005	CBR-IC3/1	CD50 (ICAM-3)	Mouse IgG1, κ
E8	353106	HA58	CD54	Mouse IgG1, κ
E9	311308	JS11	CD55	Mouse IgG1, κ
E10	398803	QA18A21	CD56 (NCAM)	Mouse IgG1, κ
E11	330905	TS2/9	CD58 (LFA-3)	Mouse IgG1, κ
E12	313906	BL-CD6	CD6	Mouse IgG1, κ
F1	336406	VI-PL2	CD61	Mouse IgG1, κ
F2	336008	HAE-1f	CD62E	Mouse IgG1, κ
F3	304806	DREG-56	CD62L	Mouse IgG1, κ
F4	304906	AK4	CD62P (P-Selectin)	Mouse IgG1, κ
F5	353004	H5C6	CD63	Mouse IgG1, κ
F6	399504	S18012C	CD64	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
B7	393003	A16085H	CD309 (VEGFR2)	Mouse IgG1, κ
B8	303106	WM59	CD31	Mouse IgG1, κ
B9	320806	1D11	CD314 (NKG2D)	Mouse IgG1, κ
B10	348406	RS38E	CD317 (BST2, Tetherin)	Mouse IgG1, κ
B11	324106	67A4	CD324 (E-Cadherin)	Mouse IgG1, κ
B12	350805	8C11	CD325 (N-Cadherin)	Mouse IgG1, κ
C1	339204	6-434	CD328 (Siglec-7)	Mouse IgG1, κ
C2	303404	WM53	CD33	Mouse IgG1, κ
C3	324306	4FR6D3	CD334 (FGFR4)	Mouse IgG1, κ
C4	331908	9E2	CD335 (NKp46)	Mouse IgG1, κ
C5	325108	P44-8	CD336 (NKp44)	Mouse IgG1, κ
C6	325208	P30-15	CD337 (NKp30)	Mouse IgG1, κ
C7	343506	581	CD34	Mouse IgG1, κ
C8	324406	24D2	CD340 (erbB2/HER-2)	Mouse IgG1, κ
C9	326606	CH3A4A7	CD344 (Frizzled-4)	Mouse IgG1, κ
C10	333406	E11	CD35	Mouse IgG1, κ
C11	314906	TREM-26	CD354 (TREM-1)	Mouse IgG1, κ
C12	353904	1D12	CD365 (Tim-1)	Mouse IgG1, κ
D1	345006	F38-2E2	CD366 (Tim-3)	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
F7	310906	FN50	CD69	Mouse IgG1, κ
F8	344004	AD2	CD73 (Ecto-5'- nucleotidase)	Mouse IgG1, κ
F9	326808	LN2	CD74	Mouse IgG1, κ
F10	341404	CB3-1	CD79b (Ilgβ)	Mouse IgG1, κ
F11	303803	QA18A37	CD8	Mouse IgG1, κ
F12	305208	2D10	CD80	Mouse IgG1, κ
G1	349506	5A6	CD81 (TAPA-1)	Mouse IgG1, κ
G2	342104	ASL-24	CD82	Mouse IgG1, κ
G3	305308	HB15e	CD83	Mouse IgG1, κ
G4	326408	17G10.2	CD85g (ILT7)	Mouse IgG1, κ
G5	333008	ZM4.1	CD85k (ILT3)	Mouse IgG1, κ
G6	336906	VIM5	CD87	Mouse IgG1, κ
G7	354104	A59	CD89	Mouse IgG1, κ
G8	301008	RPA-T8	CD8a	Mouse IgG1, κ
G9	312106	HI9a	CD9	Mouse IgG1, κ
G10	328110	5E10	CD90 (Thy1)	Mouse IgG1, κ
G11	336108	VIMD2	CD93	Mouse IgG1, κ
G12	305506	DX22	CD94	Mouse IgG1, κ
H1	305608	DX2	CD95 (Fas)	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
D2	355306	9E8	CD367 (CLEC4A)	Mouse IgG1, κ
D3	363204	m1b9	CD36L1 (SCARB1, SR-BI)	Mouse IgG1, κ
D4	397103	S17015A	CD38	Mouse IgG1, κ
D5	328208	A1	CD39	Mouse IgG1, κ
D6	334308	5C3	CD40	Mouse IgG1, κ
D7	303706	HIP8	CD41	Mouse IgG1, κ
D8	303906	HIP1	CD42b	Mouse IgG1, κ
D9	343204	CD43-10G7	CD43	Mouse IgG1, κ
D10	338808	BJ18	CD44	Mouse IgG1, κ
D11	304008	HI30	CD45	Mouse IgG1, κ
D12	323108	CC2C6	CD47	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
H2	338406	NK92.39	CD96 (TACTILE)	Mouse IgG1, κ
H3	336308	VIM3b	CD97	Mouse IgG1, κ
H4	318806	122	CD270 (HVEM, TR2)	Mouse IgG1, κ
H5	360804	22-19-12	CXCL16	Mouse IgG1, κ
H6	346404	MHD1-314	Delta-like protein 1 (DLL1)	Mouse IgG1, κ
H7	346506	MHD4-46	Delta-like protein 4 (DLL4)	Mouse IgG1, κ
H8	307106	JD3	DR3 (TRAMP)	Mouse IgG1, κ
H9	352904	AY13	EGFR	Mouse IgG1, κ
H10	357003	K152D10	GPR19	Mouse IgG1, κ
H11	358204	CG4	GPR56	Mouse IgG1, κ
H12	342604	3D12	HLA-E	Mouse IgG1, κ

Plate 3

Well	Cat #	Clone	Description	Isotype
A1			Blank	
A2	316508	MHK-49	Ig light chain κ	Mouse IgG1, κ
A3	314508	MHM-88	IgM	Mouse IgG1, κ
A4	351606	Y9A2	Integrin α9β1	Mouse IgG1, κ
A5	346904	MHJ2-523	Jagged 2	Mouse IgG1, κ
A6	346603	TDA3	Ksp37 (FGFBP2)	Mouse IgG1, κ
A7	367004	13.8	LY6G6D (C6orf23)	Mouse IgG1, κ
A8	367608	590H11G1E3	MERTK	Mouse IgG1, κ
A9	327606	W7C6	MSC	Mouse IgG1, κ
A10	330806	W4A5	MSC and NPC (W4A5)	Mouse IgG1, κ
A11	327306	W8B2	TNAP	Mouse IgG1, κ
A12	363903	TCC16	MUC-13	Mouse IgG1, κ
B1	346706	5D12	NKp80	Mouse IgG1, κ
B2	352106	MHN1-519	Notch 1	Mouse IgG1, κ
B3	345406	MHN3-21	Notch 3	Mouse IgG1, κ
B4	349004	MHN4-2	Notch 4	Mouse IgG1, κ
B5	327706	57D2	NPC (57D2)	Mouse IgG1, κ
B6	317208	NT-7	CD352 (NTB-A)	Mouse IgG1, κ

Well	Cat #	Clone	Description	Isotype
E1	362004	2E1B02	CD130 (gp130)	Mouse IgG2a, κ
E2	348506	BV9	CD144 (VE-Cadherin)	Mouse IgG2a, κ
E3	369604	BNI3	CD152 (CTLA-4)	Mouse IgG2a, κ
E4	337508	TX24	CD155 (PVR)	Mouse IgG2a, λ
E5	312606	DX27	CD158b (KIR2DL2/L3, NKAT2)	Mouse IgG2a, κ
E6	370805	SA065C3	TMEM8A	Mouse IgG2a, κ
E7	356004	K041E5	CD186 (CXCR6)	Mouse IgG2a, κ
E8	357206	K036C2	CD192 (CCR2)	Mouse IgG2a, κ
E9	353204	G043H7	CD197 (CCR7)	Mouse IgG2a, κ
E10	358904	L053E8	CD199 (CCR9)	Mouse IgG2a, κ
E11	330106	9E9A8	CD209 (DC-SIGN)	Mouse IgG2a, κ
E12	372304	W15177A	CD217	Mouse IgG2a, κ
F1	800318	3F4	CD230 (Prion)	Mouse IgG2a, κ
F2	311106	ML5	CD24	Mouse IgG2a, κ
F3	348606	UIC2	CD243 (MDR-1)	Mouse IgG2a, κ
F4	302706	BA5b	CD26	Mouse IgG2a, κ
F5	357504	19F2	CD269 (BCMA)	Mouse IgG2a, κ
F6	392305	W15145C	CD282 (TLR2)	Mouse IgG2a, κ

Well	Cat #	Clone	Description	Isotype
B7	342504	LNI-17	PSMA (FOLH1)	Mouse IgG1, κ
B8	357804	2A2	ROR1	Mouse IgG1, κ
B9	347604	5G6	Siglec-10	Mouse IgG1, κ
B10	347703	S7.7	CD328 (Siglec-7)	Mouse IgG1, κ
B11	347104	7C9	Siglec-8	Mouse IgG1, κ
B12	351504	K8	Siglec-9	Mouse IgG1, κ
C1	355203	8E11	SSEA-5	Mouse IgG1, κ
C2	327406	W5C5	SUSD2	Mouse IgG1, κ
C3	306708	IP26	α/β T Cell Receptor	Mouse IgG1, κ
C4	331210	B1	TCR γ/δ	Mouse IgG1, κ
C5	354004	9F4	Tim-4	Mouse IgG1, κ
C6	351104	MIH61	TLT-2	Mouse IgG1, κ
C7	367203	C9	TM4SF20	Mouse IgG1, κ
C8	358704	TRA-2-49/6E	TRA-2-49	Mouse IgG1, κ
C9	358803	TRA-2-54/2J	TRA-2-54	Mouse IgG1, κ
C10	322806	1B4	TSLPR (TSLP-R)	Mouse IgG1, κ
C11	356204	9D9F9	VEGFR-3 (FLT-4)	Mouse IgG1, κ
C12	400214	MOPC-173	Mouse IgG2a, κ isotype control	Mouse IgG2a, κ
D1	358604	15C4	LOX-1	Mouse IgG2a, κ

Well	Cat #	Clone	Description	Isotype
F7	312806	HTA125	CD284 (TLR4)	Mouse IgG2a, κ
F8	354704	H037G3	CD301 (CLEC10A)	Mouse IgG2a, κ
F9	354204	201A	CD303 (BDCA-2)	Mouse IgG2a, κ
F10	354504	12C2	CD304 (Neuropilin-1)	Mouse IgG2a, κ
F11	340304	509f6	CD307e (FcRL5)	Mouse IgG2a, κ
F12	356704	SHM33	CD323 (JAM3)	Mouse IgG2a, κ
G1	371204	108-17	CD357 (GITR)	Mouse IgG2a, κ
G2	336206	5-271	CD36	Mouse IgG2a, κ
G3	355404	15E2	CD369 (Dectin-1/CLEC7A)	Mouse IgG2a, κ
G4	353804	8F9	CD370 (CLEC9A/DNGR1)	Mouse IgG2a, κ
G5	353604	50C1	CD371 (CLEC12A)	Mouse IgG2a, κ
G6	304206	UCHL1	CD45RO	Mouse IgG2a, κ
G7	327910	NKI-M9	CD51	Mouse IgG2a, κ
G8	304708	p282 (H19)	CD59	Mouse IgG2a, κ
G9	343106	CD7-6B7	CD7	Mouse IgG2a, κ
G10	334106	CY1G4	CD71	Mouse IgG2a, κ
G11	326008	CD84.1.21	CD84	Mouse IgG2a, κ
G12	344304	S5/1	CD88 (C5aR)	Mouse IgG2a, κ
H1	339106	Cr24.1	CD355 (CRTAM)	Mouse IgG2a, κ

Well	Cat #	Clone	Description	Isotype
D2	375503	S18015E	CD20	Mouse IgG2a, κ
D3	318904	QA19A22	CD52	Mouse IgG2a, κ
D4	375206	S18022G	CD117 (c-kit)	Mouse IgG2a, κ
D5	394004	S16016B	CD133	Mouse IgG2a, κ
D6	367304	7.13	APCDD1 (DRAPC1)	Mouse IgG2a, κ
D7	344506	MIH26	CD272 (BTLA)	Mouse IgG2a, κ
D8	360604	L263G8	CD198 (CCR8)	Mouse IgG2a, κ
D9	358304	K097F7	CCRL2	Mouse IgG2a, κ
D10	328506	CBR-IC2/2	CD102 (ICAM-2)	Mouse IgG2a, κ
D11	327808	58XB4	CD104	Mouse IgG2a, κ
D12	355004	G077F6	CD124 (IL-4Rα)	Mouse IgG2a, κ

Well	Cat #	Clone	Description	Isotype
H2	324706	1B4C3	erbB3/HER-3	Mouse IgG2a, κ
H3	357704	K102B9	FPR3 (FPRL2)	Mouse IgG2a, κ
H4	357304	14G2a	Ganglioside GD2	Mouse IgG2a, κ
H5	363306	K07JP05	GPR83	Mouse IgG2a, κ
H6	311406	W6/32	HLA-A,B,C	Mouse IgG2a, κ
H7	307606	L243	HLA-DR	Mouse IgG2a, κ
H8	316608	MHL-38	Ig light chain λ	Mouse IgG2a, κ
H9	348204	IA6-2	IgD	Mouse IgG2a, κ
H10	337803	MHLICR2a	IL-28RA	Mouse IgG2a, κ
H11	345203	AST-3T	integrin β5	Mouse IgG2a, κ
H12	367712	SA231A2	KLRG1 (MAFA)	Mouse IgG2a, κ

Plate 4

Well	Cat #	Clone	Description	Isotype
A1			Blank	
A2	320906	6D4	MICA/MICB	Mouse IgG2a, κ
A3	327506	W3D5	SUSD2	Mouse IgG2a, κ
A4	348304	MHN2-25	Notch 2	Mouse IgG2a, κ
A5	363804	NY18	TACSTD2 (TROP2)	Mouse IgG2a, κ
A6	317307	OKT3	CD3	Mouse IgG2a, κ
A7	371306	3B2/TA8	CD99	Mouse IgG2a, κ
A8	372704	A15153G	TIGIT (VSTM3)	Mouse IgG2a, κ
A9	400314	MPC-11	Mouse IgG2b, κ isotype control	Mouse IgG2b, κ
A10	345804	hC3aRZ8	C3AR	Mouse IgG2b, κ
A11	362104	13E11	CCX-CKR (CCRL1)	Mouse IgG2b, κ
A12	371504	S-HCL-3	CD11c	Mouse IgG2b, κ
B1	310404	AH9R7	CD129 (IL-9 R)	Mouse IgG2b, κ
B2	339506	HP-MA4	CD158 (KIR2DL1/S1/S3/S5)	Mouse IgG2b, κ
B3	320608	8F1/CXCR1	CD181 (CXCR1)	Mouse IgG2b, κ
B4	310706	5E8	CD193 (CCR3)	Mouse IgG2b, κ
B5	353410	G034E3	CD196 (CCR6)	Mouse IgG2b, κ
B6	350306	51.1	CD1d	Mouse IgG2b, κ

Well	Cat #	Clone	Description	Isotype
E1	401320	MG3-35	Mouse IgG3, κ isotype control	Mouse IgG3, κ
E2	308305	CARL-1	CD255 (TWEAK)	Mouse IgG3, κ
E3	330406	MC-813-70	SSEA-4	Mouse IgG3, κ
E4	401609	MM-30	Mouse IgM, κ isotype control	Mouse IgM, κ
E5	341206	BY55	CD160	Mouse IgM, κ
E6	359612	HNK-1	CD57	Mouse IgM, κ
E7	305106	G10F5	CD66b	Mouse IgM, κ
E8	330610	TRA-1-60-R	TRA-1-60-R	Mouse IgM, κ
E9	400408	RTK2071	Rat IgG1, κ isotype control	Rat IgG1, κ
E10	347304	9-4D2-1E4	CD115 (CSF-1R)	Rat IgG1, κ
E11	351904	RCR-401	CD201 (EPCR)	Rat IgG1, κ
E12	400508	RTK2758	Rat IgG2a, κ isotype control	Rat IgG2a, κ
F1	307504	W18070C	CD161	Rat IgG2a, κ
F2	358404	3G7A02	CD120b	Rat IgG2a, κ
F3	308804	3F9	CD210 (IL-10 R)	Rat IgG2a, κ
F4	311906	1A1	CD267 (TACI)	Rat IgG2a, κ
F5	350106	BM16	CD294 (CRTH2)	Rat IgG2a, κ
F6	313612	GoH3	CD49f	Rat IgG2a, κ

Well	Cat #	Clone	Description	Isotype
B7	363504	S-HCL-1	CD22	Mouse IgG2b, κ
B8	352604	B6.220	CD220	Mouse IgG2b, κ
B9	306604	HIR2	CD235ab	Mouse IgG2b, κ
B10	318706	T5-39	CD258 (LIGHT)	Mouse IgG2b, κ
B11	329706	29E.2A3	CD274 (B7-H1, PD-L1)	Mouse IgG2b, κ
B12	331806	162.1	CD319 (CRACC)	Mouse IgG2b, κ
C1	303206	FUN-2	CD32	Mouse IgG2b, κ
C2	324206	9C4	CD326 (EpCAM)	Mouse IgG2b, κ
C3	332008	5D3	CD338 (ABCG2)	Mouse IgG2b, κ
C4	360204	9B9	CD368 (CLEC4D)	Mouse IgG2b, κ
C5	304108	HI100	CD45RA	Mouse IgG2b, κ
C6	310204	MEM-55	CD45RB	Mouse IgG2b, κ
C7	328010	NKI-SAM-1	CD49e	Mouse IgG2b, κ
C8	342304	ASL-32	CD66a/c/e	Mouse IgG2b, κ
C9	337904	24	CD85h (ILT1)	Mouse IgG2b, κ
C10	333708	GHI/75	CD85j (ILT2)	Mouse IgG2b, κ
C11	305406	IT2.2	CD86	Mouse IgG2b, κ
C12	371404	VIM15b	CD92	Mouse IgG2b, κ
D1	331104	8F11-M16	CXCR7	Mouse IgG2b, κ

Well	Cat #	Clone	Description	Isotype
F7	337704	MKT5.1	CD85a (ILT5)	Rat IgG2a, κ
F8	338706	42D1	CD85d (ILT4)	Rat IgG2a, κ
F9	410708	M1310G05	IgG Fc	Rat IgG2a, κ
F10	321204	FIB504	Integrin β7	Rat IgG2a, κ
F11	372604	S15046E	XCR1	Rat IgG2a, κ
F12	376103	W18100A	CD360 (IL-21R)	Rat IgG2a, κ
G1	337004	NC-08	Podoplanin	Rat IgG2a, λ
G2	400636	RTK4530	Rat IgG2b, κ isotype control	Rat IgG2b, κ
G3	338606	TUGh4	CD132 (common γ chain)	Rat IgG2b, κ
G4	359106	J418F1	CD195 (CCR5)	Rat IgG2b, κ
G5	357403	A161A1	CD4	Rat IgG2b, κ
G6	341604	2A9-1	CX3CR1	Rat IgG2b, κ
G7	400808	RTK2118	Rat IgM, κ isotype control	Rat IgM, κ
G8	330312	MC-631	SSEA-3	Rat IgM, κ
G9			Blank	
G10			Blank	
G11			Blank	
G12			Blank	
H1			Blank	

Well	Cat #	Clone	Description	Isotype
D2	327206	DOR7D2A4	Delta Opioid Receptor	Mouse IgG2b, κ
D3	366404	L205G1	Dopamine Receptor D1	Mouse IgG2b, κ
D4	356804	SHM16	EphA2	Mouse IgG2b, κ
D5	334610	AER-37 (CRA-1)	FcεRIα	Mouse IgG2b, κ
D6	352504	7B11	GARP (LRRC32)	Mouse IgG2b, κ
D7	330208	JM7A4	CD215 (IL-15Rα)	Mouse IgG2b, κ
D8	322008	31G4D8	Lymphotoxin beta receptor (LT-βR)	Mouse IgG2b, κ
D9	359004	K125H4	MRGX2	Mouse IgG2b, κ
D10	347504	MIH24	CD254 (TRANCE, RANKL)	Mouse IgG2b, κ
D11	324006	CUB1	CD318 (CDCP1)	Mouse IgG2b, κ
D12	621607	A17188B	CD279 (PD-1)	Mouse IgG2b, κ

Well	Cat #	Clone	Description	Isotype
H2			Blank	
H3			Blank	
H4			Blank	
H5			Blank	
H6			Blank	
H7			Blank	
H8			Blank	
H9			Blank	
H10			Blank	
H11			Blank	
H12			Blank	



LEGENDScreen™ kits are manufactured by **BioLegend**

8999 BioLegend Way
San Diego, CA 92121
United States

For contact information, as well as a list of world-wide BioLegend offices and distributors, please visit: biolegend.com/en-us/contact