


Purified anti-human FOXP3 Antibody



SeqIF™ (sequential immunofluorescence) staining on COMET™ of Purified anti-FOXP3 (clone 259D, yellow) on formalin-fixed paraffin-embedded human head and neck squamous cell carcinoma (HNSCC) tissue at 10 µg/mL. Alexa Fluor™ Plus 647 Goat anti-Mouse IgG antibody (Lunaphore, Cat. No. DR647MS) was used as a secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.