


**Brilliant Violet 421™ anti-mouse F4/80 Antibody**



C57BL/6 mouse frozen spleen section was fixed with 4% paraformaldehyde (PFA) for ten minutes at room temperature and blocked with 5% FBS plus 5% rat/mouse serum for 30 minutes at room temperature. Then the section was stained with 2 µg/ml of anti-mouse F4/80 (clone BM8) Brilliant Violet 421™ (Blue) and anti-mouse CD8a (clone 53-6.7) Alexa Fluor® 647 (red) overnight at 4°C. The image was captured with a 10X objective.