


Purified anti-human/mouse/rat CD278 (ICOS) Antibody



SeqIF[™] (sequential immunofluorescence) staining on COMET[™] of Purified anti-CD278 (ICOS) (clone C398.4A, yellow) on formalin-fixed paraffin-embedded human tonsil tissue at 5 µg/mL. Alexa Fluor[®] 647 Goat anti-Armenian hamster IgG antibody was used as secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.