


Alexa Fluor[®] 594 anti-human CD45 Antibody



SeqIF[™] (sequential immunofluorescence) staining on COMET[™] of Purified anti-HLA-DR (clone LN3, yellow) on formalin-fixed paraffin-embedded human tonsil tissue at 2.5 $\mu\text{g}/\text{mL}$. Alexa Fluor[™] Plus 647 Goat anti-Mouse IgG antibody (Lunaphore, Cat. No. DR647MS) was used as a secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.