


Purified anti- β -Amyloid, 1-11 Antibody



IHC staining of anti- β -Amyloid, 1-11 antibody (clone NAB 228) on formalin-fixed paraffin-embedded human Alzheimer's disease neocortical brain tissue. After antigen retrieval using 70% Formic Acid for 20 minutes, the tissue was incubated with the primary antibody at 10 μ g/ml for 24 hours at 4°C. BioLegend's Ultra Streptavidin (USA) HRP Detection Kit was used for detection followed by hematoxylin counterstaining.