


Purified anti-CD99 (Ewing's Sarcoma Marker) Antibody



IHC staining of anti-CD99 antibody (clone O13) on formalin-fixed paraffin embedded human Ewing sarcoma tissue. Following antigen retrieval using BioLegend's Retrieve-All Antigen Unmasking System 1, the tissue was incubated with the primary antibody for 20 minutes at room temperature at 2.7 µg/ml. BioLegend's Ultra Streptavidin (USA) HRP Kit (Multi-Species) was used for detection, followed by counterstaining with hematoxylin.