


Purified anti-Cytokeratin Antibody



IHC staining of anti-Cytokeratin antibody (clone OSCAR) on formalin-fixed paraffin-embedded human bladder tissue. Following antigen retrieval using BioLegend's Sodium Citrate Retrieval, the tissue was incubated with the primary antibody at 1 $\mu\text{g}/\text{ml}$ for 20 minutes at room temperature. BioLegend's Ultra Streptavidin (USA) HRP Detection Kit (Multi-Species) was used for detection followed by counterstaining with hematoxylin.