


Brilliant Violet 510™ anti-mouse/human CD45R/B220 Antibody



C57BL/6 mouse frozen spleen section was fixed with 4% paraformaldehyde (PFA) for ten minutes at room temperature and blocked with 5% FBS plus 5% rat/mouse serum for 30 minutes at room temperature. Then the section was stained with 1 µg/ml of anti-mouse/human CD45R/B220 (clone RA3-6B2) Brilliant Violet 510™ (green), anti-mouse CD8a (clone 53-6.7) Brilliant Violet 421™ (blue) and anti-mouse Ly-6G (clone 1A8) Alexa Fluor® 647 (red) overnight at 4°C. The image was captured with a 10X objective.