


Purified anti-Neurofilament H (NF-H) Antibody



IHC staining of anti-Neurofilament H (NF-H) antibody (clone SMI 38) on formalin-fixed paraffin-embedded rat brain tissue. After antigen retrieval using Retrieve-All 3, the tissue was incubated with the primary antibody at 5 $\mu\text{g}/\text{ml}$ for 24 hours at 4°C. BioLegend's Ultra Streptavidin (USA) HRP Detection Kit was used for detection followed by hematoxylin counterstaining.