


Brilliant Violet 421[™] anti-human CD73 Antibody



C57BL/6 mouse frozen spleen section was fixed with 4% paraformaldehyde (PFA) for ten minutes at room temperature and blocked with 5% FBS for 30 minutes at room temperature. Then the section was stained with 10 µg/ml of CD73 (clone TY/11.8) Brilliant Violet 421[™] (blue), CD3 (clone 17A2) Alexa Fluor[®] 594 (red), and B220 (clone RA3-6B2) Alexa Fluor[®] 647 (green) overnight at 4°C. The image was captured by 10X objective.