


Treg-Protector™ (anti-ARTC2 Nanobody)



Axis label: Y-% Live Treg cells.

- Cells from control mouse maintained at 4°C.
- ▨ Cells from control mouse incubated 1 hour at 37°C.
- Cells from mouse treated with Treg-Protector™, maintained at 4°C.
- ▤ Cells from mouse treated with Treg-Protector™, incubated 1 hour at 37°C.

The bar graph depicts the fraction of live Tregs recovered from mice treated with (+) or without (-) Treg-Protector™. The cells were kept at 4°C (solid bars) or at 37°C for 1 hour (hatched bars).