


Purified anti-human CD45RO Antibody



SeqIF[™] (sequential immunofluorescence) staining on COMET[™] of Purified anti-CD45RO (clone UCHL1, yellow) on formalin-fixed paraffin-embedded human tonsil tissue at 2.5 µg/mL. Alexa Fluor[™] Plus 647 Goat anti-Mouse IgG antibody (Lunaphore, Cat. No. DR647MS) was used as a secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.