


Purified anti-IDO1 Antibody



SeqIF[™] (sequential immunofluorescence) staining on COMET[™] of Purified anti-IDO1 (clone W16073A, yellow) on formalin-fixed paraffin-embedded human head and neck squamous cell carcinoma (HNSCC) tissue at 10 $\mu\text{g}/\text{mL}$. Alexa Fluor[™] Plus 647 Goat anti-Rat IgG antibody (Lunaphore, Cat. No. DR647RT) was used as secondary antibody. Nuclei were counterstained with DAPI (blue). Tissue underwent an all-in-one dewaxing and antigen retrieval preprocessing.