

Alexa Fluor® 647 anti-mouse CD20 Antibody



C57BL/6 mouse frozen spleen section was fixed with 4% paraformaldehyde (PFA) for 10 minutes at room temperature and blocked with 5% FBS plus 5% Rat Serum for 1 hour at room temperature. Then the section was stained with 10 µg/ml of CD20(clone SA275A11) Alexa Fluor®A647(red), 5µg/ml of B220(clone RA3-6B2) Alexa Fluor®A488(green) and 5 µg/ml of CD3 (clone 17A2) Alexa Fluor® 594 (blue) at 4°C overnight. The image was captured by 10X objective.